

Our Clinical Experience with the Nucleus CI532™ Electrode Array

Skřivan J., Jurovčík M., Okluský M., Bauer L., Myška P.
Department of Ear, Nose and Throat, The Second Faculty of Medicine,
University Hospital Motol, Prague, Czechia




The Nucleus CI532 electrode array was used in a clinical practice in a paediatric population.

Main features of the electrode:

1. 60% less volume compared to CI512 (Contour Advance)
2. Protection of fine cochlear structures during insertion due to the atraumatic design
3. Consistent perimodiolar position and close to spirale ganglion cells –
 - To optimise hearing outcome
 - To decrease energy consumption


On the right pictures a cross cut of the cochlea showing comparison between CI512 Contour Advance electrode (top) and CI532 Slim Modiolar electrode (bottom).
Pictures courtesy of Cochlear Corp™




Our group

Period: 3 years 9 months (6/2014 – 3/2018)
Children: 124
Gender: males 69, females 55
Age: 9 months – 16 years, average 18 months, median 24 months


Graph on the right shows a distribution of bilateral simultaneous, sequential, unilateral left / right in the whole group.
Table below shows average age of all children groups in a surgery, in case of a unilateral surgery an average age in the first and the second surgery respectively and an intersurgical interval



Type of surgery	Number of children	Average age / surgery (months)	Average age / 1 st surgery (months)	Average age / 2 nd surgery (months)	Intersurgical interval (mo)
Bilateral simultaneous	59	14,5	N/A	N/A	0
Bilateral sequential	20	N/A	28	47	19
Unilateral	45	54	N/A	N/A	N/A



Pie chart showing device distribution according to the manufacturer



Pie chart showing various types of used Nucleus™ electrodes

Complications

Due to a small diameter and flexibility of the array, the electrode has a tendency to curl inside the cochlea if a surgical technique is not correct. In our group, it occurred in two cases with the CI532 (1,5%) and in one case in a revision surgery with the CI512. Possible reasons for the tip foldover are summarized:

- Incorrect wing position during insertion
- Moving and non stable wing during insertion
- Under- or overinsertion of the polymer sheath in the basal turn
- Intracochlear changes (adhesions, obliteration)

How to detect a tip foldover intraoperatively?


- by electrophysiology (spread of excitation test SOAE, electric field imaging EFI)
- by imaging (cone-beam CT)



Postoperative high resolution CT scan, coronary plane, of a bilaterally implanted child. The electrode array in the left ear is in a good position, while a tip foldover is present in the right cochlea.

Conclusion:

The CI532 slim electrode array represents an electrode of choice in children with congenital deafness and physiological cochleas due to its perimodiolar position and atraumatic insertion



Jiri.Skrivan@fnmotol.cz